

Potentiostat

For measuring electrochemical devices and components

Cyclic voltammetry, constant voltage and constant current measurements

Impedance spectroscopy

Zahner Zennium

- 1 channel
- Equipped with high impedance probe and compliance voltage booster
- Galvanostatic output ranges
 - ± 100 nA to ± 2.5 A
 - ± 1 nA to ± 0.5 A
- Compliance voltage up to ± 120 V
- Input impedance up to 1000 T Ω
- Frequency range 10 μ Hz - 4 MHz

